


# Optical Activity

Identifying and distinguishing enantiomers is inherently difficult, since their physical and chemical properties are largely identical. Fortunately, a nearly two hundred year old discovery by the French physicist Jean-Baptiste Biot has made this task much easier. This discovery disclosed that the right- and left-handed enantiomers of a chiral compound perturb plane-polarized light in opposite ways. This perturbation is unique to chiral molecules, and has been termed **optical activity**.

Plane-polarized light is created by passing ordinary light through a polarizing device, which may be as simple as a lens taken from polarizing sun-glasses. Such devices transmit selectively only that component of a light beam having electrical and magnetic field vectors oscillating in a single plane. The plane of polarization can be determined by an instrument called a **polarimeter**, shown in the diagram below.



Monochromatic (single wavelength) light, is polarized by a fixed polarizer next to the light source. A sample cell holder is located in line with the light beam, followed by a movable polarizer (the analyzer) and an eyepiece through which the light intensity can be observed. In modern instruments an electronic light detector takes the place of the human eye. In the absence of a sample, the light intensity at the detector is at a maximum when the second (movable) polarizer is set parallel to the first polarizer ( $\alpha = 0^\circ$ ). If the analyzer is turned  $90^\circ$  to the plane of initial polarization, all the light will be blocked from reaching the detector.

András Szilágyi has created a nice animation, illustrating various kinds of polarized light. This site may be examined by [Clicking Here](#) .

Chemists use polarimeters to investigate the influence of compounds (in the sample cell) on plane polarized light. Samples composed only of achiral molecules (e.g. water or hexane), have no effect on the polarized light beam. However, if a single enantiomer is examined (all sample molecules being right-handed, or all being left-handed), the plane of polarization is rotated in either a clockwise (positive) or counter-clockwise (negative) direction, and the

analyzer must be turned an appropriate matching angle,  $\alpha$ , if full light intensity is to reach the detector. In the above illustration, the sample has rotated the polarization plane clockwise by  $+90^\circ$ , and the analyzer has been turned this amount to permit maximum light transmission.

The observed rotations ( $\alpha$ ) of enantiomers are opposite in direction. One enantiomer will rotate polarized light in a clockwise direction, termed **dextrorotatory** or (+), and its mirror-image partner in a counter-clockwise manner, termed **levorotatory** or (-). The prefixes dextro and levo come from the Latin *dexter*, meaning right, and *laevus*, for left, and are abbreviated *d* and *l* respectively. If equal quantities of each enantiomer are examined, using the same sample cell, then the magnitude of the rotations will be the same, with one being positive and the other negative. To be absolutely certain whether an observed rotation is positive or negative it is often necessary to make a second measurement using a different amount or concentration of the sample. In the above illustration, for example,  $\alpha$  might be  $-90^\circ$  or  $+270^\circ$  rather than  $+90^\circ$ . If the sample concentration is reduced by 10%, then the positive rotation would change to  $+81^\circ$  (or  $+243^\circ$ ) while the negative rotation would change to  $-81^\circ$ , and the correct  $\alpha$  would be identified unambiguously. Since it is not always possible to obtain or use samples of exactly the same size, the observed rotation is usually corrected to compensate for variations in sample quantity and cell length. Thus it is common practice to convert the observed rotation,  $\alpha$ , to a **specific rotation**,  $[\alpha]$ , by the following formula:

$$\text{Specific Rotation} = [\alpha]_D = \frac{\alpha}{l * c} \quad \begin{array}{l} \text{where } l = \text{cell length in dm, } c = \text{concentration in g/ml} \\ D \text{ is the } 589 \text{ nm light from a sodium lamp} \end{array}$$

Compounds that rotate the plane of polarized light are termed **optically active**. Each enantiomer of a stereoisomeric pair is optically active and has an equal but opposite-in-sign specific rotation. Specific rotations are useful in that they are experimentally determined constants that characterize and identify pure enantiomers. For example, the lactic acid and carvone enantiomers discussed earlier have the following specific rotations.

|   |  |
|---|--|
| Carvone from caraway: $[\alpha]_D = +62.5^\circ$          | this isomer may be referred to as (+)-carvone or <i>d</i> -carvone         |
| Carvone from spearmint: $[\alpha]_D = -62.5^\circ$        | this isomer may be referred to as (-)-carvone or <i>l</i> -carvone         |
| Lactic acid from muscle tissue: $[\alpha]_D = +2.5^\circ$ | this isomer may be referred to as (+)-lactic acid or <i>d</i> -lactic acid |
| Lactic acid from sour milk: $[\alpha]_D = -2.5^\circ$     | this isomer may be referred to as (-)-lactic acid or <i>l</i> -lactic acid |

A 50:50 mixture of enantiomers has no observable optical activity. Such mixtures are called **racemates** or racemic modifications, and are designated ( $\pm$ ). When chiral compounds are created from achiral compounds, the products are racemic unless a single enantiomer of a chiral co-reactant or catalyst is involved in the reaction. The addition of HBr to either cis- or trans-2-butene is an example of racemic product formation (the chiral center is colored red in the following equation).



Chiral organic compounds isolated from living organisms are usually optically active, indicating that one of the enantiomers predominates (often it is the only isomer present). This is a result of the action of chiral catalysts we call enzymes, and reflects the inherently chiral nature of life itself. Chiral synthetic compounds, on the other hand, are commonly racemates, unless they have been prepared from enantiomerically pure starting materials.

There are two ways in which the condition of a chiral substance may be changed:

1. A racemate may be separated into its component enantiomers. This process is called **resolution**.
2. A pure enantiomer may be transformed into its racemate. This process is called **racemization**.

Source : <http://www2.chemistry.msu.edu/faculty/reusch/VirtTxtJml/sterism2.htm#isom9>